

Sectioning with Single-View Structured Illumination

Charles A. DiMarzio
ECE, MIE, BioE
Northeastern University

April 2019

With Help From Zac Hoffman, Mahsa Azizi, Kivanc Kose
This work was supported in part by the the National Science Foundation (award number CBET-1510281).

Agenda



- Research Overview
- Modern Microscopy
- Structured Illumination (SIM)
- Random Patterns
- Sectioning with SIM
- SISIM: Single-Image SIM


Our Current Research




- Multi-Modal Microscopy
- Light and Sound
- Structured Illumination
- Collagen Orientation
- Stepwise 3-PEF in Melanin
- Lidar




Multi-Modal




Light and Sound




SIM



Collagen



Melanin



Lidar

Early Microscopes



- Compound Microscope (Jansen, 1590)
- Simple Microscope ($m=300$) (Leeuwenhoek, early 1600s)
- Physiological Observation (Hooke 1665)
- Diffraction Theory (Abbe, 1860)
- Diffraction–Limited Imaging (Spencer, mid 1880s)


Modern Microscopy



- What's so Modern?
Microscopy has been around since 1590...
- ... But a Lot Has Happened in the Last Few Decades
- Three Reasons why the Time is Right
 - Illumination Sources (From Tungsten to Lasers, LEDs)
 - Fast, Low-Cost Computers (and Cameras, *etc.*)
 - Chemistry (Molecular Tags)

What Is Sectioning?

Defocusing a camera makes an object blurry;
 Δz depends on Variation with x , ...
but contrast is an issue.




In-Focus Image



Out-Of-Focus Image


What Is Sectioning?

Defocusing a confocal microscope makes an object disappear





Judy Newmark (Warner Group), Bill Warger

Varying Spatial Frequencies



- SIM = Structured Illumination Microscopy
- Uses an Optical Fourier Transform
- Can Improve Resolution by 2X
- Can Also Provide **Sectioning**



Structured Illumination



- Image as a Product in the Field (Spatial) Plane

$$\text{Image} = \text{Illumination} \times \text{Transmission}$$

- Low Pass Filter


$$f_x < \frac{NA}{\lambda}$$

$$\text{Image} = (\text{Illumination} \times \text{Transmission}) * \text{Filter}$$


- Pupil (Spatial Frequency) Domain

$$\text{IMAGE} = (\text{ILLUMINATION} * \text{TRANSMISSION}) \times \text{FILTER}$$

How It Works




“DC” or Plane-Wave; Pupil as a Low-Pass Filter




Offset Illumination: Multiply in the Image; Convolve in the Pupil


SIM Sectioning Concept




- High frequency modulation pattern in focal plane
- Blurred pattern above and below
- Resolution dependent on frequency and NA



Fourier Domain



Layout



635nm, 18X, 0.4 NA

Random Modulation

Hoffman and DiMarzio: Structured illumination microscopy using random intensity...




Fig. 9 Wide-field *in vivo* image at the surface.




Fig. 11 Wide-field *in vivo* image at depth.




Fig. 10 CRIL *in vivo* image showing the stratum corneum.




Fig. 12 CRIL *in vivo* image showing the stratum granulosum.

SIM Sectioning Experiment


Phases = $[0^{\text{deg}}, 120^{\text{deg}}, 240^{\text{deg}}]$

$I_{\text{sec}} =$


$$\sqrt{(I_1 - I_2)^2 + (I_1 - I_3)^2 + (I_2 - I_3)^2}$$

or


$$I_{AC} = I_1 e^{i0} + I_2 e^{i2\pi/3} + I_3 e^{i4\pi/3}$$




Layer 1:
Ground Glass
Layer 2:
Gel ($n \approx 1.33$)
Layer 3:
Resolution Chart




Spiral Hilbert Transform




real(analytic signal , after spiral hilbert)




imag(analytic signal , after spiral hilbert)



abs (analytic signal , after spiral hilbert)






CCD Captures Modulated Image, $\Gamma(x, y)$.

↓
Estimate Unmodulated Image by Notch
Filtering AC Components, $R(x, y)$.


↓
Subtract Unmodulated Image from Modulated,
 $\Gamma(x, y) - R(x, y)$

↓
Demodulate Using Spiral Hilbert Transform
to Produce Sectioned Image $R_{IF}(x, y)$.



Sectioning at Depth

SISIM Phase 1 SISIM Phase 2 SISIM Phase 3




Widefield

SIM 3 Phase

SISIM Sum


Human Skin



Motion Artifacts

Human Skin at $\approx 50\mu\text{m}$

Widefield




SIM 3 Phase

SISIM

SISIM Sum
Registered


Scale Bar $\approx 10\mu\text{m}$

Registration Algorithm




- Removes modulation pattern first.
- Apply registration to just the specimen.
- Align all images before sectioning.
- Add three phases after they are aligned and sectioned
- images to produce high-quality sectioning.

Registration Results




Widefield

All Images
 $\approx 750\mu\text{m}$ Square
See Videos



SIM 3 Phase



Registered SISIM

Summary



- SIM Sectioning Artifacts Using 3 Phases
 - Refraction
 - Motion
- SISIM Uses a Single Image
- Adding Images Recovers SNR
- Registration on Image Allows Recovery of SNR

Thanks



- Funding: NSF CBET-1510281
- Graduate Students
 - Dr. Zachary Hoffman, Ph.D. 2018
 - Ms. Mahsa Azizi, Ph.D. Student
- Collaborators
 - Dr. Kivanc Kose (MSKCC)
 - Dr. Milind Rajadhyaksha (MSKCC)